

# Canfit

## TIPS and TRICKS



## Correcting a Shape Created By-Measurement

[www.vorum.com](http://www.vorum.com)

Canada/USA Locations - (800) 461-4353

Other Global Locations - (604) 321-7277

***I have entered the measurement data into the software, but the shape doesn't look right. What do I do?***

Sometimes after you have selected a library shape and have entered all of your patient's measurements into the software, the resultant shape is undesirable. Why might this happen?

1. The library shape may not suit the patient for whom you have entered the measurements.
  - √ Try a different library shape and see if the results are better.
2. While a shape was being added to the library, a landmark was placed incorrectly. This causes the measurements to be applied at an inappropriate level.
  - √ As above, try another library shape.
3. One or more of the measurements are inaccurate. This can be due to many factors, such as a stretched-out measuring tape, or a transcription error.
  - √ Try removing the measurements you entered into the software. To do this: "un-check" the white boxes, one by one (See Image 1). Only the measurements that are "checked" will be applied to the shape.




Image 1

If you need further assistance with your *Canfit*<sup>™</sup> products, please contact the Vorum support department. Email: [support@vorum.com](mailto:support@vorum.com) Phone: 1-604-321-7277